



Be Right™



## Filling Solution, Internal, Polarographic DO, 50 mL

**Product #:** 2759126  
NZD Price (Incl. GST): Contact Hach  
Not in stock - ships within 4-6 weeks

Filling solution for Polarographic Dissolved Oxygen probe, Potassium Chloride/Glycerine, 50 mL bottle

---

### Specifications

Description :	5130, 5197000
Formula:	KCl & Glycerine
Platform :	Solution
Quantity:	50 mL
Solution Type:	Refill

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)