



## Hypochlorite Test Kit, Model CN-HRDT



**Product #:** 2687100  
**NZD Price (Incl. GST):** Contact Hach  
**Not in stock - ships within 4-6 weeks**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

### Digital Titrator - Single Parameter

Utilizing the Hach Digital Titrator, this kit is used to determine hypochlorite, high range (Bleach) solution concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 5 - 15% as Cl<sub>2</sub>: 50 - 150 g/L as Cl<sub>2</sub>. Contains 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

**Designed for on-site analysis of bulk hypochlorite deliveries or bulk storage solutions at drinking water or chloramination facilities**

**Verify concentration of hypochlorite solutions at time of delivery**

**Kit contains all required reagents, solutions, and glassware**

---

### Specifications

Case Style:	G
Method:	10100
Method Name:	Digital Titrator/ Thiosulfate
Model:	CN-HRDT
Number of tests:	100
Parameter:	Hypochlorite, high range - (Bleach) As Cl <sub>2</sub>
Platform :	Digital Titrator
Range:	5 - 15 % = 50-150 g/L Cl <sub>2</sub>
Ship Wt. (lbs):	3
Smallest Increment:	0.05%

---

### What's in the box?

Includes four reagents, flask, clippers, digital titrator, delivery tube, pipet, pipet tips, instruction sheet, and carrying case.