



Be Right™



APA6000 Copper Standard 2 (high range)

Product #: 2756453
NZD Price (Incl. GST): Contact Hach
Not in stock - ships within 4-6 weeks

APA6000 Copper Standard 2 (high range), 10.0 mg/L, 1L bottle

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)